
HDU-Flow97

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SAFETY INSTRUCTIONS

Before use please read and consider the safety instructions below carefully

- Read the entire manual of the *HDM97Pocket* carefully before using the sensor.
- Keep the instrument away from unauthorised persons.
- Only use the sensor as reference device
- Never use the sensor as a replacement of sensors of a medical device.
- Do not use the sensor in conjunction with a medical device during treatment.
- Only use the instrument in a dry environment and do not touch it with wet hands.
- Make sure that no liquids intrude inside the sensor or it's connector.
- To avoid damage and misreading always confirm that the sensor is connected in the flow direction visible at the sensor.
- Make sure that the sensor filter at the inlet is clean.
- If the acquired values seem to be not believable, make sure that the sensor is not defective
- Prevent electrostatic discharge on the connectors. This can lead to substantial damage of the instrument. Make sure to be completely discharged before touching the connectors or cables connected to the instrument.
- Make sure that the instrument never heats over 60°C / 140°F. Prevent direct sunlight.
- Never sterilize the instrument using an autoclave. Danger of explosion of the battery!
- Never open the sensor. There are no parts inside you can repair.

Technical Data

The flow sensors are available in following versions. Other ranges on request.

HDU-Flow97	100.2000	50.6000	5.1750
Measuring Range H2O at 22°C / 71.6°F	100 ... 2000 ml/min	50 ... 6000 ml/min	5 ... 1750 ml/min
Accuracy From measured value	100 ... 500 ml/min ± 1.5 % otherwise ± 2.5%	+/- 1,5%	+/- 1,5%
Repeatability	Better +/- 0,5%	Better +/- 0,5%	Better +/- 0,5%
Operating Temperature	-20°C ... +80°C	-10°C ... +80 °C	-20°C ... +80°C
Operating Pressure	0 ... 5 Bar bei 22°C	-0,7... 8 Bar bei 22°C	-0,1... 5 Bar bei 22°C
Measuring Position	Variable in any direction	Variable in any direction	Variable in any direction
Filter	External	External	Internal
Fluid Connector	Dialysis Connector Male	Dialysis Connector Male	Dialysis Connector Male alternativ Tubing Connector
Picture		In preparation	
Article Number	31.0110.00	31.0111.00	31.0112.00

Intended Use

The sensors can be used in conjunction with the HDM97*Pocket* for biomedical calibration applications. In conjunction with the USB-Interface-Cable the sensor can be used in laboratory applications.

Maintenance and cleaning

Rinse the sensor with RO water after use. The inlet filter must be cleaned on a regular basis.

Calibration interval

The sensors must be calibrated once per year.


Device Connector

Sensor

Interface cable for HDM97

Interface cable for USB

Connector at HDM97 Pocket

HDM97BL

Flow

HDM97BQ with Option C (Connector Conductivity Probe)

Flow Pressure Conductivity
Temperature

HDM97BQ

Flow Pressure Conductivity
Temperature

Computer - USB

With the USB-Interface-Cable the sensors can be connected to a computer USB-Interface. For data acquisition IBPView can be used.

Kontakt Contact

IBP Medical GmbH
Ikarusallee 15
D 30179 Hannover
Germany

IBP Medical Inc
4113 W. St. Charles Ave.
Phoenix, AZ 85041
USA

Phone: +49 511 651647
Fax: +49 511 652284

1-865-686-8650
1-866-243-0187

Internet: <http://www.ibpmedical.com>
E-Mail: support@ibpmedical.com

<http://www.ibpmedical.com>
america@ibpmedical.com

Manual authors: Dipl. Ing. Werner Pfingstmann

Copyright (C) 2012 IBP Medical GmbH.
All rights reserved

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The information contained in this documentation is subject to change without notice. IBP Medical GmbH, its distributors and subsidiaries take no responsibility for any errors or omissions in this document.

All brand names mentioned in this document are property of their respective owners.